

BloodSTOP® iX

CONTROLS BLEEDING FAST

Absorbable Water Soluble Hemostat

- Etherized Oxidized Regenerated Cellulose
- Class III CE Mark approval for internal, surgical use

BloodSTOP® iX is a water-soluble, bioresorbable, animal-free hemostatic matrix which resembles traditional gauze. It has Class III CE Mark approval for internal, surgical use. It meets all international standard biocompatibility criteria for prolonged-use implant materials. Also Class II CE approved for surface use, it is the perfect choice to stop any surgical bleeding.

Not just another ORC product, BloodSTOP iX is etherized oxidized regenerated cellulose. This proprietary step makes it 100% water soluble, unlike all competitor regenerated cellulose products.

How it Works

When applied during a surgical procedure, BloodSTOP® iX quickly absorbs blood and exudates, transforms into a gel to seal the wound with a protective transparent layer, actively aids in blood coagulation, and creates an environment for wound healing. The clear gel facilitates monitoring of hemostasis.

BloodSTOP® iX Absorbable Hemostat is indicated for surgical procedures to assist in the control of capillary, venous, and small arterial hemorrhage. It can be used in neurosurgery, thoracic surgery, abdominal and gastro-intestinal surgery, plastic surgery, burn treatment, oral maxillofacial surgery, ENT surgery, obstetric and gynecologic surgery, colon-rectal surgery, circumcision, vascular surgery, ulcer wounds and other general surgeries.

Clinical experience over a ten-year period has demonstrated the efficacy of BloodSTOP iX in time-to-hemostasis, reduced blood loss and material resorption time with no toxicity or other adverse effects.

- Water soluble, natural cellulose
- Acts in under a minute
- Resorbs faster than other hemostatic wound dressings
- Designed for internal and surface wounds
- Easy to use
- Biocompatible
- Highly absorbent, hydrophilic
- No animal-derived material
- Sterile, single use packaging; several sizes
- Cost effective
- Biodegradable

Product Description	Qty per box	Product Code
BloodSTOP iX hemostat 0.5"x2" (1.3x5cm)	24/box	BS-iX-27
BloodSTOP iX hemostat 2"x2" (5x5cm)	12/box	BS-iX-14
BloodSTOP iX hemostat 2"x4" (5x10cm)	12/box	BS-iX-15
BloodSTOP iX hemostat 4"x8" (10x20cm)	12/box	BS-iX-17
BloodSTOP iX hemostat 2"x14" (5x35cm)	12/box	BS-iX-20



LifeScience PLUS, Inc.

2520-A Wyandotte Street, Mountain View, CA 94043 (or)

P.O. Box 60783, Palo Alto, CA 94306

Toll Free: 1-877-587-5433 - U.S.A. | Tel: +1-650-565-8172 - International | Fax: 1-650-336-1130

Url: www.lifescienceplus.com. For additional information, contact: sales@lifescienceplus.com